



## 24 Tube Conditioning Oven\*

### Desorption Tube and Needle Conditioning Oven for Batch Cleaning and Baking Out of Thermal Desorption Tubes and Needles.

This 24 tube conditioning oven is used for batch flow conditioning of both empty and packed glass lined SS (GLT) desorption tubes as well as for the flow conditioning of the desorption tube needles. The system is designed much like the 6 port conditioning oven but contains a much larger heater block with 24 ports for tubes and needles. Four port manifold handles permit the flow conditioning of 4 tubes or needles from one ball rotameter.

Each of the six adjustable flow bubble rotameters has a flow range between 0 and 50 ml/min, one for each of the manifolds with tubes or needles to be conditioned. Each of the flows to the manifolds can be independently controlled via a microneedle valve and the flow turned off to that port when not needed. A single electrically operated solenoid valve (Gas Valve) turns off or on the gas to all the ports via the Gas Switch on the front panel. On the top of the Conditioning Oven six quick disconnects connect to 1/8" flexible nylon tubing to provide gas flow through the manifolds with desorption tubes and needles during conditioning. When the quick disconnect is removed, the gas flow is automatically closed to that port. By providing gas flow from a carrier gas such as high purity helium or nitrogen through the desorption tubes and needles while conditioning, it can be assured that no oxygen enters the GLT desorption tube which could destroy the adsorbent material. Impurities from the inside of the desorption tubes and needles are flushed through. A single 1/8" fitting on the back of the Conditioning System is provided for the attachment of the carrier gas from the source.



Conditions 24 Desorption Tubes or Needles Simultaneously

Electrical power to the conditioning system oven is provided from a standard 110 volt, 10 amp outlet. An external high temperature reset switch is located on the rear of the Oven so that if temperatures exceed 350°C on the heater block, the external reset will open and heat to the block will cease. This switch can be reset by simply pushing in the button once the block has cooled sufficiently. Two slo-blow fuses are mounted on the rear of the Conditioning System, one for the main power and one for the heater circuit. In addition a three inch fan enclosed on the back side of the Conditioning System provides a steady flow of air through the system to maintain the temperature of the enclosed electronic components to an acceptable level. If the case temperature exceeds 60°C, an internal switch will cut power to the heater circuit.

Also located on the back of the Conditioning System is a cooling rack. After the heated conditioning of the desorption tubes and needles is complete, the desorption tubes and needles are removed from the oven and placed in the cooling rack. As soon as the tubes are touchable ( 5 to 10 minutes), they are immediately capped on both ends with the stainless steel caps with PTFE seals that are available for the desorption tubes.

#### Specifications

##### Electrical Specifications

Power Requirements  
Voltage - 110 VAC, grounded  
Current - 10 amp maximum  
Fuse Sizes - 10 amp each

##### Temperature Controller

Heater Circuit, Accuracy  $\pm 0.1\%$  of full scale reading  
Temperature Range - up to 350°C  
Input - Platinum Resistance Thermometer  
Digital Readout - 3 or 4 digit  
Programmable - Autocalibration temperature controller  
Temperature protection circuitry for maximum temperature and maximum case temperature

##### Gas / Flow Specifications

6 adjustable Flow bubble rotameters  
Maximum pressure - 60 PSI  
Quick disconnects for each gas line  
Gases - either high purity helium or nitrogen

##### Weight and Dimensions

Weight - 20 pounds  
Size - 9" wide x 9" deep x 13" high

Warranty - 90 days parts and labor

Part No.	Description	Price
781056	24 Port Conditioning Oven, 6 rotameters, programmable temperature controller, 6 - 4 port desorption tube manifold handles, 1 - 4 port needle manifold handle with graphite seals	
781057	4 Port Manifold Desorption Tube Handle	
781058	4 Port Manifold Needle	
781015	Graphite Seal with Metal Insert	
781007	S. S. Cap with .040" Hole	

\* Recommended for use with the AutoDesorb System