

The Right Parts Make All the Difference

**Quick Guide for Agilent GC/MS
Replacement Parts and Supplies**

For Optimum Performance . . . use Authentic Agilent GC/MS Replacement Parts and Supplies


Want the most value from your Agilent gas chromatograph/mass spectrometer investment? Your Agilent GC/MS system will run at its very best when you choose Agilent-specified consumables and replacement parts.



Agilent's *GC/MS Resource Guide* (publication 5988-3960EN) tells how to get peak performance from GC/MS systems.

Sample Vials

- Vials, clear wide opening crimp top, write-on spot (100/pk), **5182-0543**
- Crimp caps, silver-aluminum, PTFE/silicone septum (100/pk) **5182-0552**
- 100 μ L conical glass microvolume vial inserts, pulled point (100/pk) **5183-2085**



MSD Supplies

- Ion source cleaning kit **5181-8863**
- PFTBA sample, certified (10 g) **8500-0656**
- El filament, **05972-60053**
- Electron multiplier replacement horn **05971-80103**
- Diffusion pump oil (128 mL) **6040-0819**
- Rough pump oil (1 liter) **6040-0834**



Autosampler Syringes

- 5 μ L, 23-gauge fixed needle ALS syringe (6/pk), **5182-0875**
- 10 μ L, 23-gauge fixed needle ALS syringe (6/pk), **9301-0725**



Gas Purifier

- High capacity carrier gas purifier, 1/8" fittings, removes O₂, H₂O and hydrocarbons from helium, **RMSH-2**
- Split vent kit, **RDT-1020**


Inlet Supplies

- Gold plated inlet seal, **18740-20885**
- Bleed and Temperature Optimized 11 mm septa (50/pk), **5183-4757**
- High-pressure Merlin Microseal starter kit, **5182-3442**
- 4 mm ID straight tube split liner (5/pk), **210-3003-5**
- 4 mm ID single-taper splitless liner (5/pk), **5183-4695**



Column Fittings and Ferrules

- Inlet column hex nut (2/pk), **5181-8830**
- Graphite/VespeI *inlet* ferrules for 250 μ m and smaller ID columns (10/pk) **5181-3323**
- Graphite/VespeI *inlet* ferrules for 320 μ m ID columns (10/pk), **5062-3514**
- MSD interface column nut, **05988-20066**
- Graphite/VespeI *MSD* ferrules for 250 μ m and smaller ID columns (10/pk) **5062-3508**
- Graphite/VespeI *MSD* ferrules for 320 μ m ID columns (10/pk), **5062-3506**



Learn More

www.agilent.com/chem

800 227 9770

or contact your local
Agilent office or authorized
Agilent distributor.

www.agilent.com/chem
800 227 9770

For More Information

To find out more about Agilent GC/MS supplies, our Customer Contact Center is a convenient, single point of contact for everything from product information and quotations to ordering and technical assistance. Visit online at www.agilent.com/chem or call 1-800-227-9770 in the U.S. and Canada for quick, knowledgeable help. Or contact your local Agilent office or authorized Agilent distributor.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2002
All Rights Reserved

Printed in USA March 4, 2002
5988-5428ENUS



Agilent Technologies