

Manufacturer	Instrument	Manufacturer No.	Burle Model No	S.I.S. Part No.	Price
Sciex	ELAN 5000, 500, 250, 6000, 6000A, 6100	N810-1069	4816B	30111	
	ELAN5000A	WE02-0328	4880	30146	
	API III LCMS	014032A	4879	30145	
	API 100, 150, 200, 300, 3000, 4000, 5000, Q-TRAP	020048B	4822B	30118	
	4000 Q-TRAP	-	4822C	33440	
	API2000, API3200	-	4899	35119	
Shimadzu	QP1000, QP2000	-	4870E	30139	
	QP5000, QP5050A, QP8000	-	4797	34225	
Stanford Research	RGA100, 200, 300	-	4220	33929	
Teledyne	3DQ	-	5902 MAGNUM	32336	
Thermo	1015 -On-axis, 3000	30004-60160	4750G	32102	
	1015 -Off axis	30004-60160	4751G	30081	
	3000, 3100, 3200, 4000,4500, 4600	01504-60110	4751G	30081	
	TSQ@70 with 20kV dynode, TSQ@ 700	94011-60052	4715	30071	
	TSQ70, ITMS, INCOS@50, INCOS@XL	40007-98020	4752G	30082	
	Sola ICP	-	4870V	30140	
	PPINICI 4000, 4500, 4600, 5100 (Base)	40007-98020	4752G	30082	
	OWA @ (before 8/81), 1020, 5100	30004-60160	4750G	32102	
	OWA (after 8/81), 1020, 5100 (Base)	30004-60050	4752G	30082	
	CH7, CH8	422097	4733	33866	
	ITD™ 700, 800, ITS40™, Magnum™	-	4765G	30088	
			5905 MAGNUM	34002	
	GCQ, LCQ, SSQ7000, TSQ7000, PolarisQ, TSQ Quantum	-	5903 MAGNUM™	32906	
	GCQ, LCQ, SSQ7000, TSQ7000, PolarisQ, TSQ Quantum	-	5903 MAGNUM™ (Entire assembly)	33119	
	LCQ Decca	-	5905 MAGNUM	34002	
Thermo Jerrel Ash	POEMS	-	4870V	30140	
UTI	Analog RGA (unmounted)	100C	4717	30074	
	Analog RGA (mounted)	100C	4776	30102	
	Pulse Counting RGA (mounted)	100C	4876	30144	
V&F	CIMS 500	-	4823G	30121	
Varian	Auger Analyzer	-	4720	30075	
	1200	-	5904 MAGNUM™	33328	
	Saturn (Pre-2000), Saturn I, II, III (mounted)	-	4765G	30088	
	Saturn 2000	-	4765GM	32754	
	Saturn 2100, 2200	-	4715	30071	
	Saturn™ I & II (unmounted), 4000 GC/MS	-	4765G	30088	
Vestec	201 (off-axis)	-	4782	30105	
	201 (on-axis) +ion mode only	-	4783	33871	
VG/	Replacement for Philips X919VL	-	7010	30150M	
Thermo VG Scientific	Replacement for Philips X919BL	-	7010M	30151	
	Replacement for Philips B310VL	-	7012	33871	
	Replacement for Philips X812VL	-	7013C/B1	30155	
	QX200, SQ300, SQ400	-	4771	30095	
	SX200, SX300, SQ300	MGR4073B	4774	30099	
	ZAB1F, 7070	-	4773G	30098	
VG Analytical	ZAB, 7070	-	4773G	30098	
VG Elemental	Genesis, PQ-3, PQ-EXCEL	2027814	4870V	30140	
VG Microtech	CLAM2 Hemispherical Analyzer (single channel detection)	-	7010	30150M	
	CLAM2 Hemispherical Analyzer (triple channel detection)	-	4800B1 (4x19mm)	30106	
VG Quadrupoles/	Quad Mass Soec (analog)	-	4774	30099	
Thermo ONIX	VGQ, VG Smart IQ	-	4769	30090	
VSW	Mass Analyst Quadrupole Analog Type	-	4771	30095	

PPINICI™ is a trademark of Thermo Corp. The Thermo 3200 PRINICI multiplier has the Z bracket mounted on the bottom plate reversed. Be sure to add the suffix -3200 to all PPINICI multipliers for the 3200 and S.I.S. will reverse the Z bracket on the multiplier base.


The Thermo OWA® used 2 styles of multipliers. OWA® is a registered trademark of Thermo Corp.

\* Adapters - several multipliers require a special adaptor to retrofit existing instruments. These adapters are reusable and need only be purchased the first time a Burle multiplier is installed.

\*\* The cost of conversion is covered by the return of your old multiplier. See the Thermo multiplier page of this catalog for more details.

† These Channeltrons contain a grid or screen over the channel entrance.

Channeltron® is a registered trademark of Burle Electro Optics Corp.



CHANNELTRON® 4752

When "Thermo" is indicated, it is intended to imply "Thermo Fisher Scientific Corp." \*When "Burle" is indicated, it is intended to imply "Photonis USA, Inc."