

# SGE Ferrules for Perkin Elmer and Shimadzu GC D49

## Perkin Elmer

| Material  | Column I.D.                  | Size of Nut | Ferrule I.D. | Pack size | Part No. | Price |
|---|------------------------------|-------------|--------------|-----------|----------|-------|
| <b>15% Graphite / 85% Vespel Ferrules</b>   |                              |             |              |           |          |       |
| <b>For Injectors &amp; Detectors at atmospheric pressure e.g. FID</b>             |                              |             |              |           |          |       |
|   | 0.1 - 0.25mm                 | 1/16"       | 0.4mm        | 10        | 072663   |       |
|   | 0.1 - 0.25mm                 | 1/8"        | 0.4mm        | 10        | 0726703  |       |
|   | 0.32mm                       | 1/16"       | 0.5mm        | 10        | 072654   |       |
|   | 0.32mm                       | 1/8"        | 0.5mm        | 10        | 0726702  |       |
|   | 0.45 - 0.53mm                | 1/16"       | 0.8mm        | 10        | 072655   |       |
|   | 0.45 - 0.53mm                | 1/8"        | 0.8mm        | 10        | 072671   |       |
|   | for 1/8" O.D. packed columns | 1/8"        | 1/8"         | 10        | 072669   |       |
|   | for 1/4" O.D. packed columns | 1/4"        | 1/4"         | 10        | 072667   |       |
| <b>100% Graphite Ferrules</b>   |                              |             |              |           |          |       |
| <b>Injectors &amp; Detectors at atmospheric pressure e.g. FID (Not for GC/MS)</b> |                              |             |              |           |          |       |
|   | 0.1 - 0.32mm                 | 1/16"       | 0.5mm        | 10        | 072627   |       |
|   | 0.1 - 0.32mm                 | 1/8"        | 0.5mm        | 10        | 072624   |       |
|   | 0.45 - 0.53mm                | 1/16"       | 0.8mm        | 10        | 072626   |       |
|   | 0.45 - 0.53mm                | 1/8"        | 0.8mm        | 10        | 0726280  |       |
|   | for 1/8" O.D. packed columns | 1/8"        | 1/8"         | 10        | 072622   |       |
|   | for 1/4" O.D. packed columns | 1/4"        | 1/4"         | 10        | 072621   |       |
| <b>SilTite™ Metal Ferrules</b>  |                              |             |              |           |          |       |
| <b>For GC/MS Interface Connection (Starter kit)</b>                               |                              |             |              |           |          |       |
|   | 0.1 - 0.25mm                 | -           | 0.4mm        | 10*       | 073200   |       |
|   | 0.32mm                       | -           | 0.5mm        | 10*       | 073201   |       |
|   | 0.53mm                       | -           | 0.8mm        | 10*       | 073202   |       |
| <b>Replacement SilTite™ Metal Ferrules</b>  |                              |             |              |           |          |       |
| <b>For All Injectors and Detectors - Including GC/MS</b>                          |                              |             |              |           |          |       |
|   | 0.1 - 0.25mm                 | -           | 0.4mm        | 10        | 073220   |       |
|   | 0.32mm                       | -           | 0.5mm        | 10        | 073221   |       |
|   | 0.53mm                       | -           | 0.8mm        | 10        | 073222   |       |
| <b>Replacement SilTite™ Nuts</b>  |                              |             |              |           |          |       |
|   | -                            | -           | -            | 5         | 073224   |       |

\*Includes 10 ferrules, 2 SilTite™ nuts.

## Shimadzu

| Material  | Column I.D.             | Ferrule I.D. | Pack size | Part No. | Price |
|---|-------------------------|--------------|-----------|----------|-------|
| <b>15% Graphite / 85% Vespel Ferrules</b>   |                         |              |           |          |       |
| <b>For QP5000-I GC/MS System</b>  |                         |              |           |          |       |
|   | 0.1 - 0.25mm            | 0.4mm        | 10        | 0726563  |       |
|   | 0.32mm                  | 0.5mm        | 10        | 0726564  |       |
| <b>For QP5000-II &amp; QP5050 GC/MS System</b>  |                         |              |           |          |       |
|   | 0.1 - 0.25mm            | 0.4mm        | 10        | 0726561  |       |
|   | 0.32mm                  | 0.5mm        | 10        | 0726562  |       |
| <b>For 2010 Injector &amp; GC/MS Connection</b>                                       |                         |              |           |          |       |
|   | 0.1 - 0.25mm            | 0.4mm        | 10        | 072663   |       |
|   | 0.32mm                  | 0.5mm        | 10        | 072654   |       |
|   | 0.53mm                  | 0.8mm        | 10        | 072655   |       |
| <b>100% Graphite Ferrules</b>   |                         |              |           |          |       |
| <b>For Injectors &amp; Detectors at atmospheric pressure e.g. FID (Not for GC/MS)</b> |                         |              |           |          |       |
|   | 0.1 - 0.32mm            | 0.5mm        | 10        | 0726080  |       |
|   | 0.45 - 0.53mm           | 0.8mm        | 10        | 0726082  |       |
|   | 5mm O.D. Packed Columns | 5.0mm        | 10        | 0726001  |       |
| <b>For 2010 Injector</b>  |                         |              |           |          |       |
|   | 0.1 - 0.32mm            | 0.5mm        | 10        | 072635   |       |
|   | 0.45 - 0.53mm           | 0.8mm        | 10        | 072636   |       |
| <b>SilTite™ Metal Ferrules</b>  |                         |              |           |          |       |
| <b>For GC/MS Interface Connection - Direct Interface (Starter kit)</b>                |                         |              |           |          |       |
|   | 0.1 - 0.25mm            | 0.4mm        | 10*       | 073204   |       |
| <b>For GC/MS Interface Connection - Jet Separator &amp; 2010 (Starter kit),</b>       |                         |              |           |          |       |
|   | 0.1 - 0.25mm            | 0.4mm        | 10*       | 073200   |       |
|   | 0.32mm                  | 0.5mm        | 10*       | 073201   |       |
|   | 0.53mm                  | 0.8mm        | 10*       | 073202   |       |
| <b>For All Injectors, Jet Separator (Starter kit), except 2010</b>                    |                         |              |           |          |       |
|   | 0.1 - 0.25mm            | 0.4mm        | 10        | 073350   |       |
|   | 0.32mm                  | 0.5mm        | 10        | 073351   |       |
|   | 0.53mm                  | 0.8mm        | 10        | 073352   |       |
| <b>Replacement SilTite™ Metal Ferrules</b>  |                         |              |           |          |       |
| <b>For GC/MS Interface Connection</b>   |                         |              |           |          |       |
|   | 0.1 - 0.25mm            | 0.4mm        | 10        | 073220   |       |
|   | 0.32mm                  | 0.5mm        | 10        | 073221   |       |
|   | 0.53mm                  | 0.8mm        | 10        | 073222   |       |
| <b>For GC/MS Interface Connection - Direct Interface (For Part No. 073204)</b>        |                         |              |           |          |       |
|   | 0.1 - 0.25mm            | 0.4mm        | 10        | 073227   |       |
| <b>Replacement SilTite™ Nuts</b>  |                         |              |           |          |       |
| <b>For 2010 GC/MS System</b>  |                         |              |           |          |       |
|   | -                       | -            | 5         | 073224   |       |

\*Includes 10 ferrules, 2 SilTite™ nuts.