

D50 SGE GC Ferrules - Instrument Cross Reference

Agilent

Description	Part No.
Column Connection - 15% Graphite-85% Vespel ferrules for all injectors and detectors. Not for GC/MS	
For 0.10 - 0.25 mm id Columns*	073109
For 0.32 mm id Columns*	073111
For 0.45 - 0.53 mm id Columns*	073113
For 1/8" OD Packed Columns	072669
For 1/4" OD Packed Columns	072667

Column Connection - 100% Graphite ferrules for all injectors and detectors. Not for GC/MS

For 0.10 - 0.32 mm id Columns*	072635
For 0.45 - 0.53 mm id Columns*	072636
For 1/8" OD Packed Columns	072602
For 1/4" OD Packed Columns	072601

Column Connection - 15% Graphite-85% Vespel ferrules. For GC/MS

For 0.10 - 0.25 mm id Columns	072663
For 0.32 mm id Columns	072654
For 0.45 - 0.53 mm id Columns	072655

*Ferrule body style (short) fits standard Agilent injector/detector capillary nut.

Perkin Elmer

Description	Part No.
Column Connection - 15% Graphite-85% Vespel ferrules	
For 1/16" nut, 0.10 - 0.25 mm id Columns	072663
For 1/16" nut, 0.32 mm id Columns*	072654
For 1/16" nut, 0.45 - 0.53 mm id Columns	072655
For 1/8" nut, 0.22 - 0.25 mm id Columns	0726703
For 1/8" nut, 0.32 mm id Columns	0726702
For 1/8" OD Packed Columns	072669
For 1/4" OD Packed Columns	072667

Column Connection - 100% Graphite ferrules

For 1/16" nut, 0.10 - 0.32 mm id Columns	072627
For 1/16" nut, 0.45 - 0.53 mm id Columns	072626
For 1/4" nut, 1/4" OD Packed Columns	072621
For 1/8" nut, 0.22 - 0.32 mm id Columns	072624
For 1/8" nut, 0.45 - 0.53 mm id Columns	0726280
For 1/8" nut, 1/8" OD Packed Columns	072622

Shimadzu

Description	Part No.
Column Connection - 15% Graphite-85% Vespel ferrules	
0.10 - 0.25 mm id Columns	0726080
0.45 - 0.53 mm id Columns	0726082
For 5 mm Packed Columns	0726001

Column Connection - 100% Graphite ferrules

QP5000-I for 0.10 - 0.32 mm id Columns	0726563
QP5000-I for 0.32 mm id Columns	0726564
QP5000-II & QP5050 for 0.22 - 0.32 mm id Columns	0726561
QP5000-II & QP5050 for 0.32 mm id Columns	0726562

*When "Thermo" is indicated, it is intended to imply "Thermo Fisher Scientific Corp."

Thermo*

Description	Part No.
Column Connection - 15% Graphite - 85% Vespel ferrules - not suitable for 8000 series GC	
For 0.10 - 0.25 mm id Columns*	073109
For 0.32 mm id Columns*	073111
For 0.45 - 0.53 mm id Columns*	073113
For 1/8" OD Packed Columns	072669
For 1/4" OD Packed Columns	072667

Column Connection - 100% Graphite ferrules for all injectors and detectors. Not for GC/MS

For 0.10 - 0.32 mm id Columns*	072635
For 0.45 - 0.53 mm id Columns*	072636
For 1/8" OD Packed Columns	072602
For 1/4" OD Packed Columns	072601

Low Cost Capillary Conversion Nut Kit and Ferrules for Thermo Finnigan - CE Instruments/Fisons 8000

Supplied with range of ferrules for 0.15 - 0.53mm id Columns	103436
--	--------

*Ferrule body style (short) fits standard Agilent injector/detector capillary nut.

Varian

Description	Part No.
Column Connection - 15% Graphite - 85% Vespel ferrules	
For 0.10 - 0.25 mm id Columns	072663
For 0.32 mm id Columns	072654
For 0.45 - 0.53 mm id Columns	072655
For 1/8" OD Packed Columns	072669
For 1/4" OD Packed Columns	072667

Column Connection - 100% Graphite ferrules

For 0.22 - 0.32 mm id Columns	072627
For 0.45 - 0.53 mm id Columns	072626
For 1/8" OD Packed Columns	072622
For 1/4" OD Packed Columns	072621

All ferrules come in packages of 10