


E32 SIS™ Solvent Delivery Accessories

SIS Mobile Phase Filter/Sparging Frits



Part Number	Porosity	Tubing size	Max Flow Rate		Price ea.
			Liquid	Helium	
E100-2	2µm	1/16"	8 ml/min	--	
E100-10	10µm	1/16"	20 ml/min	400 ml/min	
E101-2	2µm	1/8"	8 ml/min	--	
E101-10	10µm	1/8"	20 ml/min	400 ml/min	

► *Tech Tip* For Helium sparging, use 10µm porosity.

SIS Economy Mobile Phase Filter

Use every drop of your precious mobile phase. These economically priced filters are easy to use and clean. 316 stainless steel filter in an inert, durable PTFE housing. Use with 1/8" tubing.




Part Number	Porosity	Tubing size	Filter Diameter	Price ea.
E110-2	2µm	1/8"	19mm	
E110-10	10µm	1/8"	19mm	
E111-2	2µm	1/8"	25mm	
E111-10	10µm	1/8"	25mm	

SIS Biocompatible Mobile Phase Filter

Eliminate metals from your mobile phase completely! Great for LC/MS of proteins and IC applications. Use with 1/8" O.D. tubing.

Part Number	Description	Price
E120	Biocompatible Mobile Phase Filter includes 5µ porous polyethylene filter, PEEK fittings, Tefzel ferrule	ea.
E121	Replacement filter element	Pk/2



SIS Mobile Phase Cap

Fits all standard bottles with GL-45 thread. The cap has three holes with 1/4"-28 threads for sparging, delivery, or return lines. All wetted parts are PTFE.



Part Number	Description	Price
E130	Cap with PEEK plug	
E131	Cap Kit with: 1 Cap, 3m 1/8" PTFE tubing, 1 Mobile Phase Filter 2µm, 1 Sparging Frit 10µm, 1 p/n E546 PEEK Plug	

► *Tech Tip* Use P/N E231 Flangeless Fitting and P/N E904 Tefzel ferrule for leak-free connections.