


Biphasic PicoFrit Columns

Biphasic PicoFrit® columns combine ProteoPep™ II C18 chemistry with strong cation exchange.




Part No.	Column ID (µm)	Tip Size (µm)	Bed Length (cm)	Pieces per Pkg.	Price
PFC7515-CXPP2-10	75	15	10	1	
PFC7515-CXPP2-10-3	75	15	10	3	
PFC5010-CXPP2-10	50	10	10	1	
PFC5010-CXPP2-10-3	50	10	10	3	

Self-Pack PicoFrits

Ready for loading - Perfect for Proteomics

PicoFrit® columns are ready for you to pack with your preferred polymer or silica media. Column packing is made easy by the low back-pressure frit.

- Spray directly from column to inlet
- Minimize post-column loss



Self-Pack PicoFrits® Ordering Information

Part Number	Tubing Size (OD/ID µm)	Nominal Tip ID (µm)	Frit Length (µm)	Tubing Length (cm)	Pieces per Pkg.	Price
Standard Coating						
PF360-75-30-CE-5	360/75	30 ± 2.0	350	50	5	
PF360-75-15-CE-5	360/75	15 ± 1.5	250	50	5	
PF360-75-10-CE-5	360/75	10 ± 1.0	200	50	5	
PF360-50-30-CE-5	360/50	30 ± 2.0	350	50	5	
PF360-50-15-CE-5	360/50	15 ± 1.5	250	50	5	
PF360-50-10-CE-5	360/50	10 ± 1.0	200	50	5	
No Coating						
PF360-75-30-N-5	360/75	30 ± 2.0	350	50	5	
PF360-75-15-N-5	360/75	15 ± 1.5	250	50	5	
PF360-75-10-N-5	360/75	8 ± 1.0	200	50	5	
PF360-50-30-N-5	360/50	30 ± 2.0	350	50	5	
PF360-50-15-N-5	360/50	15 ± 1.5	250	50	5	
PF360-50-10-N-5	360/50	8 ± 1.0	200	50	5	
Self-Pack PicoFrit® Assortment Kits - Includes: (1) PF360-75-8, (2) PF360-75-15, and (2) PF360-75-30						
FRITKIT-360-CE	Standard Coating					
FRITKIT-360-N	No Coating					