

| Vial Size (mm) | 8 x 40 | 12 x 32 | 12 x 32 | 12 x 32 | 12 x 32 | 12 x 32 | 12 x 32 | 12 x 32 | 15 x 45 | 15 x 45 |
|---|------------|------------|--------------|------------|---------|------------------|--------------------|---------------------|---------------------|------------|
| Vial Type | Shell Vial | Shell Vial | Crimp Finish | Snap/Crimp | ABC | 9mm Screw Thread | 8-425 Screw Thread | 10-425 Screw Thread | 13-425 Screw Thread | Shell Vial |
| Mfr. & Model | | | | | | | | | | |
| Dani | | | | | | | | | | |
| ALS 86.80 | - | - | • | • | • | • | • | - | - | - |
| ALS 39.80 | - | - | • | • | • | • | • | - | - | - |
| Dionex/Gynkotek | | | | | | | | | | |
| Gina | - | - | • | • | - | - | • | - | - | - |
| ALI-100 | - | - | - | • | - | - | • | - | - | - |
| AS 509 | - | - | - | • | - | - | • | - | - | - |
| Dynatech/Precision | | | | | | | | | | |
| WPS-3000 | - | - | • | • | • | • | • | - | - | - |
| 42 Place Tray | - | - | • | • | • | • | - | - | - | - |
| 60 Place Tray | - | - | - | - | - | - | - | - | - | - |
| LC2000 | - | - | - | - | • | • | - | - | - | - |
| 231 XL, 232 XL | - | - | • | • | • | • | • | - | - | - |
| ESA/LC Packing | | | | | | | | | | |
| Model 542 | - | - | • | • | • | • | • | - | - | - |
| Midas | - | - | • | • | • | • | • | - | - | - |
| Gilson | | | | | | | | | | |
| 233 Sample Changer | - | - | - | - | - | - | - | - | - | - |
| 231 XL/232XL/ 233XL | - | - | • | • | • | • | • | - | - | - |
| GX-271/281/ | - | - | • | • | • | • | • | - | - | - |
| Aspec XL | - | - | • | • | • | • | • | - | - | - |
| Hitachi | | | | | | | | | | |
| L2200 | - | - | • | • | • | • | • | - | - | - |
| L7250 | - | - | • | • | • | • | • | • | • | - |
| L8800 | - | - | - | • | • | • | - | - | - | - |
| L7200 | • | - | - | - | - | - | - | - | • | • |
| AS-2000 | - | - | • | • | • | • | • | - | - | - |
| AS-4000 | - | - | • | • | • | • | • | - | - | - |
| Jasco | | | | | | | | | | |
| AS-2050/2055/2057 | - | - | • | • | • | • | • | - | - | - |
| AS-2059 | - | - | • | • | • | • | • | - | - | - |
| AS-1555/1555-10 | - | - | • | • | • | • | • | - | - | - |
| AS-1559 | - | - | • | • | • | • | • | - | - | - |
| 851-AS, AS-900 | - | - | • | • | • | • | • | - | - | - |
| LC 800/900 Series | - | - | - | - | - | - | • | - | - | - |
| Leap Technologies | | | | | | | | | | |
| CTC DI PAL | - | - | • | • | • | • | • | - | - | - |
| CTC Combi PAL | - | - | • | • | • | • | • | - | - | - |
| CTC A200E | - | - | • | • | • | • | • | - | - | - |
| CTC HTS PAL | - | - | • | • | • | • | • | - | • | - |
| CTC HTS Twin PAL | - | - | • | • | • | • | • | - | - | - |
| CTC HTC PAL | - | - | • | • | • | • | • | - | - | - |
| CTC LC mini PAL | - | - | • | • | • | • | • | - | - | - |
| CTC A2000S | - | - | • | • | • | • | • | - | - | - |
| A200LC | - | - | • | • | • | • | - | - | - | - |
| Perkin Elmer | | | | | | | | | | |
| Claris 500, 600 | - | - | • | • | • | • | • | - | - | - |
| Autosystem | - | - | • | • | • | • | • | - | - | - |
| Autosystem XL | - | - | • | • | • | • | • | - | - | - |
| Series 200 | - | - | • | • | • | • | • | - | - | - |
| Integral 4000 | - | - | • | • | - | - | - | - | - | - |
| ISS 100, ISS 200 | - | - | • | • | • | • | • | - | - | - |
| LC 600 42 Place Tray | - | - | • | • | • | • | - | - | - | - |
| Phillips/ Pye UNICAM | | | | | | | | | | |
| 4710 | - | - | • | • | - | - | • | - | - | - |
| 4700LC-GC/S4/S8 | - | - | - | - | - | - | - | - | - | - |
| 4247 | - | - | • | • | - | - | - | - | - | - |
| C-XP | - | - | - | - | - | - | - | - | • | • |
| Shimadzu GC | | | | | | | | | | |
| AOC-20/AOC-20s | - | - | - | - | - | - | • | • | • | • |
| AOC-8B | - | - | - | - | - | - | • | • | - | - |
| AOC-5000 | • | - | - | - | - | - | - | - | - | - |
| AOC-14/AOC-17 | - | - | - | - | - | - | • | • | • | • |
| Shimadzu HPLC | | | | | | | | | | |
| SIL-HTa/SIL-HTc | • | - | - | • | • | • | • | • | • | - |
| SIL-10ADVP | • | - | - | • | • | • | • | • | • | - |
| SIL-10A/SIL-10Ai/SIL-10Ap | • | - | - | • | • | • | • | • | • | - |
| SIL-6B/SIL-9A/SIL-8A | - | - | - | - | - | - | • | • | • | - |
| SIL-10AD | - | - | - | - | - | - | • | • | • | - |
| Promis | • | - | • | • | • | • | • | - | - | - |
| Spark Holland | | | | | | | | | | |
| Reliance | - | - | • | • | • | • | • | - | - | - |
| Midas | - | - | • | • | • | • | • | - | - | - |
| Marathon | • | - | • | • | • | • | • | - | - | - |
| Triathon | • | - | • | • | • | • | • | - | - | - |
| Thermo Scientific (division of Thermo Fisher Scientific) | | | | | | | | | | |
| AS1000/AS3000/AS3500 | - | - | • | • | - | - | • | - | - | - |
| TriPlus AS/HS/DUO | - | - | • | • | - | - | • | - | - | - |
| Surveyer Plus / Plus Lite | - | - | • | • | - | - | • | - | - | - |
| Accela | - | - | • | • | - | - | • | - | - | - |
| 8875/8880 | - | - | • | • | - | - | • | - | - | - |
| Varian GC | | | | | | | | | | |
| CP-8400/CP-8410 | - | - | • | • | • | • | • | - | - | - |
| 8200CX | - | - | - | - | - | - | - | - | - | - |
| Ultra GC/MS | - | - | - | - | • | • | - | - | - | - |
| 8100/8200 | - | - | • | • | • | • | • | • | - | - |
| Varian LC | | | | | | | | | | |
| 9100/9090/9095 | - | - | • | • | - | - | • | - | - | - |
| ProStar 400/410 | - | - | • | • | - | - | - | - | - | - |
| ProStar 420/430 | - | - | • | • | - | - | - | - | - | - |
| Marathon | • | - | • | • | • | • | • | - | - | - |
| Waters LC | | | | | | | | | | |
| Acquity | - | - | - | • | • | • | • | • | - | - |
| Alliance | • | - | - | • | • | • | • | - | • | • |
| Breeze | • | - | - | • | • | • | • | - | - | - |
| 717Plus | • | - | - | • | • | • | • | - | • | • |
| 2700 Sample Manager | - | - | - | • | • | • | • | - | - | - |
| 710/712/715 | - | - | • | - | - | - | - | - | • | • |