

GCQ® Filament Repair

S.I.S. can repair the Thermo MAT GCQ filament (Thermo #70001-60287). Due to the complexity of the design, this filament is more difficult to repair and therefore the number of times it can be repaired may be limited.



Part No.	Description	Price ea.
FILREP	Repair of Filament	
	6 or more	

DEP Probe Tips

The DEP Probe Tips for the GCQ (Thermo #119322-0001) are manufactured by S.I.S. In addition, our shops can repair your old filaments. We can also custom manufacture these probe tips to your specific requirements.



Part No.	Description	Price ea.
FINDEPGCQ	GCQ-DEP Probe Tip (Quantity of 6 or more)	
FILREP	Repair of Probe Tips (Quantity of 6 or more)	

Fuses

Part No.	Thermo #	Description	Price Pkg of 5
218160	00006-01700	0.16A, 250V, S.B. 5x 20mm	
218250	00006-11204	0.25A, 250V, S.B., 5x20mm	
2180315	00006-04550	0.315A, 250V, S.B., 5 x 20mm	
218400	00006-05080	0.4A, 250V, S.B., 5 x 20mm	
218500	00006-11105	0.5A, 250V, S.B., 5 x 20mm	
218800	00006-06075	0.8A, 250V, S.B., 5 x 20mm	
218001	00006-11101	1.0A, 250V, S.B., 5 x 20mm	
2180125	00006-07650	1.25A, 250V, S.B., 5 x 20mm	
218016	00006-08650	1.6A, 250V, S.B., 5 x 20mm	
218025	00006-09510	2.5A, 250V, S.B., 5 x 20mm	
218025	00006-11202	2.5A, 250V, S.B., 5 x 20mm	
2183150	00006-10510	3.15A, 250V, S.B., 5 x 20mm	
217002	00006-08806	2.0A, 250V, Fast Acting	
217001	00006-07610	1.0A, 250V, Fast Acting	
217500	00006-07608	0.5A, 250V, Fast Acting	

Vespel™ Transfer Line Tip

The Vespel Transfer Line tip (96000-20026) is now manufactured at S.I.S.




Part No.	Description	Price ea.
VA800	GCQ Transfer Line Tip	

Viton® O-Rings

Part No.	Thermo #	Thickness	I.D.	Price each	Price Pkg of 10
V126	00107-10400	3/32"	1.37"	-	
V120	00107-10401	3/32"	1.0"	-	
V154	00107-10058			-	
V214	00107-11500	1/8"	1.0"	-	
K1091050 (to 550°F)	00107-10059				
K1094079 (to 600°F)					

Thermo GCQ High Temperature Direct Insertion Probe

- Maximum temperature 800° C
- Fast ramp rates
- For use on Thermo GCQ
- Direct heating of samples
- Strong, easy to use flared vials



High Temperature Probe Mounted On Probe Inlet System

Cal Gas Probe

A new Cal Gas Probe can also be used with the GCQ Probe Inlet System. See Section A for details.

Part No.	Description	Price
PP400	Probe Inlet for GCQ	
FP16	High Temperature Probe for GCQ	
PC2	High Temperature Probe Controller	
S17A	Flared Sample Vials, Package of 100	
S17B	Flared Sample Vials, Package of 1000	
S17C	Flared Sample Vials, Package of 500	
CAL3	Cal Gas Probe for GCQ	

A complete line of multipliers are also available for the GCQ from SIS. See Section A for details

*When "Thermo" is indicated, it is intended to imply "Thermo Fisher Scientific Corp"